



MSR Gold™ MiniFastFit (DJ/Club)

MSR Gold™ 300/2 MiniFastFit 1CT/4

Thanks to the specially designed lamp base and lamp holder, the MSR Gold™ MiniFastFit takes just seconds to replace, to help ensure non-stop clubbing entertainment. In addition, the innovative gold plated pins ensure excellent current transfer to prolong lamp and lamp holder lifetimes and thereby reduce lamp replacement costs. Easy handling and small compact fixture designs are allowed by the lamps compact size, while the very short arc provides a high perceived brightness and high beam intensity.

Warnings and Safety

- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.
- Lamp contains mercury.
- Manage in Accord with Disposal Laws.
- See: www.lamprecycle.org or 1-800-555-0050

Product data

General information		Initial lumen (Nom)		21300 lm
Base	PGJX28 [PGJX28]	Chromaticity Coordinate X (Nom)	288	
Operating Position	Universal [Universal]	Chromaticity Coordinate Y (Nom)	306	
Main Application	Entertainment	Correlated Color Temperature (Nom)	9300 K	
Life to 50% Failures (Nom)	750 h	Luminous Efficacy (Rated) (Min)	64 lm/W	
System Description	MiniFastFit	Luminous Efficacy (rated) (Nom)	71 lm/W	
Light technical		Color Rendering Index (Nom)	80	
Color Code	2	Operating and electrical		
Initial lumen (Min)	19170 lm	Power (Rated) (Nom)	300 W	

MSR Gold™ MiniFastFit (DJ/Club)

Lamp Current (Nom)	3.75 A
Ignition Supply Voltage (Min)	207 V

Controls and dimming

Dimmable	Yes
----------	-----

Mechanical and housing

Cap-Base Information	PGJX28
----------------------	--------

Luminaire design requirements

Bulb Temperature (Max)	1742 °F
Pinch Temperature (Max)	932 °F

Product data

Order product name	MSR Gold™ 300/2 MiniFastFit 1CT/4
EAN/UPC - Product	8718291221111
Order code	287177
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	4
Material Nr. (12NC)	928177105115
Copy Net Weight (Piece)	0.137 lb

Dimensional drawing



Product	D (max)	O	L (min)	L (max)	L	C (max)
MSR Gold™ 300/2 MiniFastFit 1CT/4	15/16 in	1/4 in	2-5/8 in	2-11/16 in	2-5/8 in	5 in

MSR Gold™ 300/2 MiniFastFit 1CT/4

