



MSD Platinum (Touring/Stage)

MSD Platinum 20 RB UNP/24

The MSD Platinum lamps open up new levels of creative freedom in entertainment lighting, with powerful, compact and brilliantly intense light. Their compact, lightweight design provides the freedom to create smaller, lighter luminaires – for use at any location on stage. And their short arc, high color temperature and innovative reflector produce a sparkling, high beam intensity for outstanding color reproduction. They also offer long and reliable lifetime, high efficacy and fast replacement times. It all adds up to an intensely exciting creative experience.

Warnings and Safety

- A lamp breaking is extremely unlikely to have any impact on your health. If a lamp breaks, ventilate the room for 30 minutes and remove the parts, preferably with gloves. Put them in a sealed plastic bag and take it to your local waste facilities for recycling. Do not use a vacuum cleaner.
- Lamp contains mercury. Manage in Accord with Disposal Laws. See: www.lamprecycle.org or 1-800-555-0050

Product data

General Information	
Base	- [-]
Operating Position	Universal [Universal]
Main Application	Stage/Touring
Life to 50% Failures (Nom)	1500 h
System Description	Reflector E21.9 TCM
Light Technical	
Color Code	2
Initial lumen (Nom)	22400 lm
Chromaticity Coordinate X (Max)	314
Chromaticity Coordinate X (Min)	280

Chromaticity Coordinate X (Nom)	296
Chromaticity Coordinate Y (Max)	316
Chromaticity Coordinate Y (Min)	290
Chromaticity Coordinate Y (Nom)	298
Correlated Color Temperature (Nom)	8000 K
Luminous Efficacy (rated) (Nom)	47 lm/W
Operating and Electrical	
Power (Rated) (Nom)	470 W
Voltage (Max)	110 V
Voltage (Min)	85 V

MSD Platinum (Touring/Stage)

Controls and Dimming

Dimmable	Yes
----------	-----

Mechanical and Housing

Cap-Base Information	na [-]
----------------------	---------

Luminaire Design Requirements

Bulb Temperature (Max)	900 °C
Pinch Temperature (Max)	350 °C

Product Data

Order product name	MSD Platinum 20 RB UNP/24
EAN/UPC - Product	8718696684221
Order code	323733
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	24
Material Nr. (12NC)	928495205308
Net Weight (Piece)	136.600 g

Dimensional drawing

D

Product	D (max)	O	C (max)
MSD Platinum 20 RB UNP/24	58.5 mm	1.2 mm	72 mm

MSD Platinum 20 RB

