



# EvoKit 2x2 HE Gen5

## EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5

The Philips EvoKit LED Retrofit Kits are an energy efficient, easy to install solution to upgrade your fluorescent troffers to LED. Compatible with both standard and narrow T-grids, they offer a simple retrofit that will improve the look of your ceiling with its architectural styling without the need to actually break the ceiling plenum. The units also come standard with dimming capabilities, making them perfect for applications such as offices, classrooms, healthcare facilities, retail space and more. The Philips EvoKit LED Retrofit Kits offer the latest advances in LED technology, resulting in quality lighting with extremely high efficacies of up to 162 lm/w. Both the 2'x2' and the 2'x4' EvoKit are manufactured with quality components and finishes, meaning a consistent, balanced lighting scheme when using both configurations in the same space. The main diffuser and slanted troffer help reduce glare and create a pleasant, uniform throw of the light. Combine the aesthetics and quality with the ease of installation, and this product can literally transform your space in minutes!

### Product data

General Information		Conduit Connection	
Lighting Technology	LED	7/8 inch knockout	
Mounting	Recessed	Mechanical and Housing	
Number Of Pieces	1	Shape	Square
		Fixture Size	2 ft x 2 ft
		Enclosure	Metal and Plastic
Light Technical		Product Data	
Luminous Flux	3,100 lumen	Order product name	EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5
Correlated Color Temperature (Nom)	5000 K	Order code	EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5
Luminous Efficacy (rated) (Nom)	155 lm/W	Numerator - Quantity Per Pack	1
Color rendering index (CRI)	80	Numerator - Packs per outer box	1
Operating and Electrical		Material Nr. (12NC)	929002715913
Power Consumption	20 W	Catalog Number	EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5
Voltages	120/277 V		

## EvoKit 2x2 HE Gen5

<b>Catalog Number Description</b>	EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5	<b>EAN/UPC - Product/Case</b>	781087165966
<b>Full product name</b>	EvoKit HE 2x2 31L 20W 850 UNV 0-10 G5		
<b>EAN/UPC - Case</b>	50781087165961		

